
CLEANING, DISINFECTION AND STERILIZATION GUIDELINES

en	4
es	10
de	16
fr	22
it	28
pt	34
cs	40
ro	46
lv	52
lt	58
da	64
hu	70
sv	76
et	82
fi	88
el	94
no	100

CLEANING, DISINFECTION AND STERILIZATION GUIDELINES

1. INTRODUCTION


The aim of this document is to establish the cleaning (manual and automatic), disinfection (manual and automatic) and sterilization procedures for devices manufactured by Biotechnology Institute S.L. (hereinafter BTI).

This applies for first-time use of multiple-use components after removal of the protective transport packaging and after subsequent uses, as well as for single-use devices that are delivered nonsterile.

With these cleaning, disinfection and sterilization guidelines, BTI provides a validated procedure to ensure clean, disinfect and sterile products. According to EN ISO 17664, it remains the responsibility of the processor to ensure that the reprocessing using cleaning and disinfectant agents, materials and equipment in the reprocessing facilities achieves the desired result as actually performed. Likewise, any deviation by the processor from the instructions provided should be properly evaluated for effectiveness and potential adverse effects.

2. PROCEDURE

The preparation for reuse of medical devices or for initial use of new devices requires a sequence of processing steps shown in this workflow:

Protecting clothes, instruments, and equipment

Whichever method is used, staff should always use suitable protective clothing (water impermeable protection gown, gloves, face/eye shield, etc.), instruments and equipment.

Cleaning agents and disinfectants

Only cleaning and disinfectant agents that meet the requirements of the applicable medical device standards and regulations in each market should be used for the cleaning and disinfection of medical devices.

BTI does not recommend any specific cleaning and/or disinfection agents but, since not all cleaning and disinfectant agents may be available worldwide, criteria for the selection are also provided.

A list of cleaning and disinfectant agents and conditions for use that BTI has used in the validation of these cleaning, disinfection and sterilization guidelines is provided.

Table 1 Cleaning and disinfectant agents and conditions for use

Method	Detergent/disinfectant (criteria for selection)	Detergent/disinfectant (supplier)	Conditions for use (Solution, Temperature, Time)
Pre-treatment	Detergent for the reprocessing of thermostable and thermolabile instruments: Alkaline detergent with enzymes and surfactants	Neodisher MediClean Forte (Dr. Weiger)	0.5% (5ml/l) / 20°C (room temperature) / 30 min
Automated cleaning	Detergent for the reprocessing of thermostable and thermolabile instruments: Alkaline detergent with enzymes and surfactants	Neodisher MediClean Forte (Dr. Weiger)	0.3% (3ml/l) / 55°C / 5 min
Manual cleaning	Detergent for the reprocessing of thermostable and thermolabile instruments: Alkaline detergent with enzymes and surfactants	Neodisher MediClean Forte (Dr. Weiger)	0,5% (5ml/l) / 40°C / 30 min
Automated disinfection	N/A	N/A	Thermal disinfection. See Table 2.
Manual disinfection	High-level disinfectant for reprocessing semi-critical heat-sensitive and reusable medical devices: ortho-Phthalaldehyde	Cidex OPA Solution (Johnson & Johnson)	0,55% (undiluted) / 20°C (room temperature) / 5 min

Materials and Equipment

Only materials and equipment (washer disinfector, autoclave, packaging, etc.) that meet the requirements of the applicable standards and regulations in each market should be used for the cleaning, disinfection and sterilization of medical devices.

BTI does not recommend any specific material or equipment but, since not all may be available worldwide, criteria for the selection are also provided.

A list of materials and equipment and conditions for use that BTI has used in the validation of these cleaning, disinfection and sterilization guidelines is provided.

Table 2 Materials and equipment and conditions for use

Materials / equipment	Criteria for selection	Materials / equipment (supplier)	Conditions for use (Temperature, Time, etc.)
Ultrasonic bath	CE marking according to Directive 98/79/ EC for IVD	Selecta HD22L 39-40kHz	40°C / 30 min / > 35 kHz
Washer-disinfector	Complies with the standard ISO 15883	Miele Unit G 7836 CD	Program parameters / Temperature / Duration Prewashing I / 10°C / 10 minutes Washing / 55°C / 5 minutes Neutralization / 10°C / 2 minutes Rinsing II / 10°C / 1 minute Thermic disinfection / 93°C / 5 minutes Drying / 110°C / 25 minutes
Packaging	Complies with the standard ISO 11607-1 or EN 868-5	Sealable pouch of cellulose porous material and plastic film according to EN 868-5	self-sealed

Pre-vacuum Autoclave	Complies with the standard ISO 17665-1 or EN ISO 13060 or EN 285	Autester STDRYPVII	132°C / 4 minutes (US method) 134°C / 5 minutes (EU method)
Containers, syringes, wipes, measuring jugs, soft bristle brushes, etc.	N/A	N/A	N/A

Cleaning and disinfectant agents and conditions for use used by BTI in the validation is described in table 1. Materials and equipment and conditions for use used by BTI in the validation is described in table 2.

2.1 PROCESSING AT THE POINT OF USE AND TRANSPORT

After each surgical procedure, all reusable devices should be cleaned of gross soiling with absorbent wipes and placed in a bath containing a cleaning solution suitable for this according to the detergent manufacturer’s instructions to prevent blood or body fluids from drying before decontamination. Devices should be kept as short as possible before cleaning and disinfection according to the detergent manufacturer’s instructions.

Where disassembly of a medical device is necessary for the cleaning and disinfection process to be effective, it shall be disassembled in accordance with the instructions for use of the medical device itself.

Single-use devices do not require this pre-treatment prior to cleaning and disinfection as long as there is no risk of cross-contamination with biological residues (e.g., accidental contact with used devices). In case of doubt, pre-treatment should be carried out as described for reusable devices.

Single-use devices or worn tools must be properly disposed after use in accordance with the general hygiene directives and legal regulations governing the proper disposal of infectious material.

If devices need to be transported to other facilities for processing, put the devices in trays and protect them with a damp cloth wipe. This will prevent any mechanical damage from occurring and biological debris from drying out.

2.2 CLEANING AND DISINFECTION

2.2.1 AUTOMATIC CLEANING AND DISINFECTION

-  Be sure that detergent is suitable for use with medical devices and washer-disinfector machine.
-  Be sure that medical devices are suitable for thermal disinfection.
-  Tools must not have contact among them to prevent possible physical damage.

1. Remove gross soiling by brushing every device or parts of devices with a soft bristle nylon brush or absorbent wipes and rinse with tap water at room temperature for at least 3 minutes in order to eliminate most of the impurities. Special attention should be paid on areas where access is difficult such as seals, angles, blind/narrow holes, etc. as well as rough surface or laser engraving.
2. Disassemble the medical device, if applicable, according to the instructions for use of the device. Clamps and other devices with joints shall be opened to avoid deposition on these parts and torque tools shall be in non-active position.
3. Place devices inside the washer – disinfector machine according to the manufacturer’s instructions. Load devices such that hinges are open and cannulations and holes can drain.
4. Add the cleaning solution suitable for this according to the detergent manufacturer’s instructions and run the cleaning and disinfection program according to the instructions of the washer-disinfector supplier.
5. After that, run the drying program and if needed, dry devices manually using filtered medical compressed air by blowing off/blowing through ensuring that every tool is completely dry.
6. Carry out a visual inspection of devices paying special attention on areas where access is difficult such as seals, angles, blind holes, etc. If traces of impurities are observed, repeat all cleaning and disinfection steps from point 1 to 5.

Cleaning and disinfectant agents and conditions for use used by BTI in the validation is described in table 1. Materials and equipment and conditions for use used by BTI in the validation is described in table 2.

2.2.2 MANUAL CLEANING

-  Be sure that detergent is suitable for use with medical devices and ultrasonic bath.
-  Tools must not have contact among them to prevent possible physical damage.
-  Do not exceed a temperature over 45°C due to the risk of protein coagulation.
-  Cleaning devices by ultrasonic bath could not be suitable for items that contain carbon of tungsten in their composition.
-  Cleaning devices by ultrasonic bath is not recommended for cleaning of rubber or plastic items as these tend to absorb ultrasounds.
-  Cleaning devices by ultrasonic bath is not recommended for cleaning items that contain an opaque layer suitable for scanning as ultrasounds can damage the layer.
-  The cleaning solution should be replaced as soon as turbidity or suspended particles being detected due to cleaning activity could be inhibited.

1. Remove gross soiling by brushing every device or parts of devices with a soft bristle nylon brush or absorbent wipes and rinse with tap water at room temperature for at least 3 minutes in order to eliminate most of the impurities. Special attention should be paid on areas where access is difficult such as seals, angles, blind/narrow holes, etc. as well as rough surface or laser engraving.
2. Disassemble the medical device, if applicable, according to the instructions for use of the device. Clamps and other devices with joints shall be opened to avoid deposition on these parts and torque tools shall be in non-active position.
3. Fully immerse the devices in an ultrasonic bath with purified water and add the cleaning solution suitable for this according to the detergent manufacturer’s instructions.

4. Maintain the devices in the ultrasonic bath according to the instructions of the supplier of the cleaning agent and the ultrasonic bath.
5. Remove devices from the ultrasonic bath and brush them again with a soft bristle nylon brush or absorbent wipes and rinse with tap water at room temperature for at least 1 minute. Special attention should be paid on areas where access is difficult such as seals, angles, blind/narrow holes, etc. as well as rough surface or laser engraving.
6. Flush all surfaces and lumens 5 times with at least 10 ml each of purified water at room temperature using a single-use syringe.
7. Carry out a visual inspection of devices paying special attention on areas where access is difficult such as seals, angles, blind holes, etc. If traces of impurities are observed, repeat all cleaning steps from point 1 to 6.

Cleaning and disinfectant agents and conditions for use used by BTI in the validation is described in table 1. Materials and equipment and conditions for use used by BTI in the validation is described in table 2.

2.2.3 MANUAL DISINFECTION

-  Be sure that disinfectant is suitable for use with medical devices.
-  Tools must not have contact among them to prevent possible physical damage.

1. Place the devices that have just been cleaned in a bath containing a disinfectant solution suitable for this according to the disinfectant manufacturer’s instructions. Be sure that devices are fully immersed and there is no contact among them.
2. Remove devices from the disinfectant solution and rinse thoroughly by immersing them completely in a large volume of potable water (e.g. 8L) for at least 1 minute.
3. Flush all surfaces and lumens 5 times with at least 20 ml each of purified water at room temperature using a single-use syringe. Repeat this process two more times to get three independent rinses.
4. After that, dry the devices manually by blowing with filtered medical compressed air and/or clean, lint-free single-use wipes as necessary, ensuring that all tools are completely dry.

- Carry out a visual inspection of devices paying special attention on areas where access is difficult such as seals, angles, blind holes, etc. If traces of impurities are observed, repeat all disinfection steps from point 1 to 5.

! If sterilization process could not be possible or recommended after disinfection, be sure to use sterilized water for rinsing in order to avoid recontamination. A final rinse using a 70% isopropyl alcohol solution can be used to speed the drying process.

! The disinfectant solution should be replaced according to the instructions for use of the supplier.

Cleaning and disinfectant agents and conditions for use used by BTI in the validation is described in table 1. Materials and equipment and conditions for use used by BTI in the validation is described in table 2.

2.3 INSPECTION, MAINTENANCE, ASSEMBLING AND PACKAGING

- Before sterilization, all medical devices should be inspected, assembled and packed in an appropriate packaging.
- As a general rule, a visual inspection without magnification will be sufficient, but good lighting conditions (e.g. 500 lux) are necessary. For this inspection at least the following should be checked: is the device visibly clean? is there any rust or corrosion? are there any broken or deformed parts? are the cutting edges dull or worn out?
Additionally, it shall be verified that the device has not reached the end of its shelf life or the recommended number of uses.
- After inspection, the devices shall be assembled and lubricated with a water-soluble lubricant, if necessary, according to the device instructions. In addition, the device shall be checked to ensure that it is correctly assembled and that it functions correctly according to the instructions for use.

! Do not use lubricants based on silicone or oils due to inactivate sterility process.

- After assembly and verification of the correct functioning of the device, it shall be placed in sterilisation pouches or blisters. Ensure that the pack is large enough to contain the device without stressing the packaging.

! The correct application of suitable wrapping is critical for the success of the sterilization.

! Be sure that devices with locks or articulations and torque tools are in non-active position, otherwise there is a possibility of breakage due to expansion caused by heat and inadequate drying after sterilisation.

! Instruments with sharps or cutting edges like drills must be packed separately in pouches or sterilised using a sterilization box.

- If a sterilisation box (e.g. autoclavable BTI surgical box) is used, position the items so that they are fixed and place the box in a sterilisation pouch of sufficient size without stressing the packaging.

Materials and equipment and conditions for use used by BTI in the validation is described in table 2.

2.4 STERILIZATION AND DRYING

- Place the sterilisation pouches or blisters in a pre-vacuum autoclave in an upright or flat position, side by side, with the transparent film surfaces facing each other. Make sure not to exceed the maximum autoclave load specified by the supplier and not to stack the pouches.
- Sterilise devices according to the following parameters:
 - 4 minutes at 132°C (US method)
 - 5 minutes at 134°C (EU method)
- After completing the sterilization cycle, let dry and cool bags inside the autoclave for at least 20 minutes, before removing and storing them.

2.5 STORAGE

Place the sealed bags in a closed and dry place. Make sure there is enough space between the instruments and the floor, walls and ceiling for circulate enough air.

Follow the instructions of the manufacturer of the packaging regarding storage conditions and expiration date of sterilized goods.

3. SYMBOLS

For a description of the symbols on the product labels and in these instructions, please refer to the MA-087 guide.

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