
# ENDORET®

## PLASMATERM H-PLUS

MANUAL OF USE AND MAINTENANCE

CE  
0123


Rx only

**WARNING:**  
Before using the PLASMATERM H-PLUS, please read this manual carefully.

	<b>en</b>	PLASMATERM H-PLUS OVEN.....	(Page 2)	
	<b>es</b>	HORNO PLASMATERM H-PLUS .....	(Página 8)	
	<b>de</b>	PLASMATERM H-PLUS OFEN.....	(Seite 14)	
	<b>fr</b>	FOUR PLASMATERM H-PLUS.....	(Page 20)	
	<b>it</b>	FORNO PLASMATERM H-PLUS .....	(Página 26)	
	<b>pt</b>	FORNO PLASMATERM H-PLUS .....	(Pagina 32)	

# ENDORET®

---

## PLASMATERM H-PLUS

### 1. INTENDED USE:

The **PLASMATERM H-PLUS** is an oven designed to heat and maintain at 37°C or 56°C Plasma Rich in Growth Factors (PRGF), within the scope of Endoret® technology. Since its use implies handling blood fractions, the unit is intended to be used exclusively by qualified health personnel in hospital or clinical environments.

### 2. DESCRIPTION:

The **PLASMATERM H-PLUS** is a robust and easy to operate unit, thanks to its internal metal casing and its exterior plastic cover which makes it highly resistant to external damage.

The oven includes a door for accessing the internal cavity. The door is easily blocked open to facilitate the manipulation of the recipients.

Inside the cavity there is a tray for placing the plasma containing recipients. The tray, made of anodised aluminium, is resistant, light, extractable and easy to clean.

The front information panel of the oven has 4 light-emitting diodes (leds) indicating at all times the operational process being used.

The working temperature inside the oven is achieved by air circulation. Hot air generated by a resistance located inside an aluminum radiator is circulated by a fan to the interior of the oven and then back to the radiator. The overall control of the equipment is carried out by a micro-processor, as well as the temperature control and protection systems.

 **WARNING:**  
Use only the tray supplied with the equipment.

#### ACCESSORIES

The unit is supplied without accessories.

 **WARNING:**  
Be sure to use recipients biocompatible and resis-

*tant to 59 °C. The maximum amount of plasma to be placed inside the oven is 100ml. Exceeding the maximum amount will not guarantee the correct heating of the plasma at the working temperature.*

#### SPARE PARTS:

In order to supply any original spare parts, the code or description of the component and the serial number of the unit should be indicated.

 **WARNING:**  
BTI shall not be held responsible for damage that may be caused by non-original spare parts. Use only original spare parts.

Picture 1

### 3. CONTRAINDICATIONS AND ADVERSE EFFECTS:

*No contraindications or adverse effects have been identified.*

*Pregnant and nursing women: No extra precautions have been identified.*

## 4. INSTALLATION

This chapter contains all the instructions needed to install the unit correctly in its working place.

### INSPECTION AND UNPACKAGING

Inspect the overall condition of the packaging, if you observe any deterioration (wet or broken) contact your distributor. After opening the packaging, extract the power cable and the unit with the tray, and verify that the list of contents attached coincides with the material delivered. If any material is missing, is not as specified or is defective, contact your distributor.

**WARNING:**

To transport the unit, pick it up by the lower part of the casing. Holding the unit by other part could damage some components of the unit (see Picture 1).

**DANGER:**

When moving or unpackaging the unit, pay attention to the tray as it comes as an independent part, not fixed to the rest of the unit. Incorrect handling when transporting could damage the unit and the user.

### LOCATION

Place the unit in a horizontal surface to avoid plasma spills. Leave a safety distance all around the equipment to ensure a safety use as well as a good air circulation. The table on which the unit is placed must have access to a power socket, at a maximum distance of 1 m. This power supply must have the characteristics required for the unit to work correctly. See section "8. Technical characteristics".

### POWER SUPPLY

The power connector is located on the rear (at the right side) including the master switch and the fuse box (see Picture 2). Before connecting the power cable to the unit, check the power supply values of your provider, (voltage, frequency, etc.) are correct. Compare these values with those assigned for the unit on the characteristics label (see Picture 3) next to the power connector or in section "8. Technical Characteristics" of this manual. The power supply must be earthed and the unit must be connected to this.

Picture 2

CHARACTERISTICS LABEL FUSE BOX SWITCH

### CHARACTERISTICS LABEL

Picture 3

### INSTALLATION

1. Once the unit has been unpacked and placed at its working place, you must ensure the safety distances are suitable according to section 4.2 and that the unit's rubber supports are all placed on the working surface.
2. Check the master switch, close to the power connector, is in the "0" position (off).
3. Connect the unit making sure that the power supply is as specified on the characteristics label or in section 8.

The unit is now installed and ready for use.

**WARNING**

Check the unit's power supply requirements, so that it is not connected to an unsuitable source.

## 5. OPERATION

To start using the unit, first you must familiarise yourself with it and with its front panel, which will give you the necessary information at all times.

### FRONT PANEL

Picture 4 shows the front panel. For a better understanding, a table is attached with all the parts that comprise it and a small explanation for each of them.

Picture 4

- 1 STATUS button
- 2 STATUS led
- 3 PROG button
- 4 LED 1
- 5 LED 2
- 6 LED 3

<b>1. Turning ON the unit</b>	Connect the unit to the power supply and press the master switch, at the rear of the unit, to "1" position (ON). The STATUS led starts blinking. Press the STATUS button in the front panel. The STATUS turns ON. Additionally, the 3 other LEDs (1, 2, 3) start blinking and then turn OFF.
<b>2. Program 1: selecting a working temperature of 37 °C</b>	Press the PROG button to select Program 1. LED1 will blink until the temperature of 37°C is reached. Once reached, LED1 will stop blinking and remain ON permanently.
<b>3. Program 2: selecting a working temperature of 56 °C</b>	Press the PROG button to select Program 2. LED2 will blink until the temperature of 56°C is reached. Once reached, LED2 will stop blinking and remain ON permanently.
<b>4. Program 3: changing from 37 °C to 56 °C for 1h.</b>	This program lasts an hour and allows changing from an initial temperature of 37 °C to a final temperature of 56 °C. Program 3 can only be accessed from Program 1. When LED1 is ON permanently, then press PROG button to select Program 3. LED3 starts blinking increasing its frequency as it approaches to the end of the hour. When the unit reaches the working temperature (56 °C) LED 2 turns ON. At the end of the hour, alarm will sound to indicate the time is off. Press the PROG button again to stop the alarm and turn off the equipment.

**Note:** If you detect any problem related to the unit functioning contact BTI Technical Service.

**FUNCTIONING OF THE UNIT**

Turn on the unit as follows:

1. Press the master switch, on the power connector, to "1" position (ON).
2. Press the STATUS button.
3. Select the desired program (see section 5.1)

After using the equipment, turn off as follows:

1. Press the STATUS button; the unit is OFF. The STATUS led turns OFF.
2. Press the master switch, on the power connector, to "0" position (OFF).

**WARNING**  
 To work with the 37°C program after using the 56°C program, remember to open the door of the oven to help lowering the temperature.  
 During the temperature transition, LED1 will be blinking and once the unit is stable at 37°C LED1 will be ON permanently.

**USE OF BIOLOGICAL MATERIAL**

The unit is designed to heat and maintain Plasma Rich in Growth Factors at 37°C or 56°C, within the scope of the Endoret® technology.

Handling of blood fractions implies the risk of presence of pathogenic microorganisms, so only professional health-care staff must use the oven and must have a suitable level of protection to ensure their safety against any contamination.

The unit is to be used only for its intended purpose and by professional healthcare staff in a hospital or clinical environment. Any use beyond the specifications of this manual is under the full responsibility of the end user.

**DANGER**  
 Take the necessary precautions and use appropriate decontaminating products to protect against any contamination from accidental spills. An important guide of biosafety can be obtained in the Laboratory Biosafety Manual, edited by the World Health Organisation. It is also necessary to comply with the current regulations related to exposure to biological agents. Check section 6. "Maintenance" of this manual.

## 6. MAINTENANCE

It is responsibility of the user to carry out the following procedures when necessary, and keep the unit in good working conditions to ensure the correct functioning of the equipment.

If repair or maintenance is needed for the unit you have acquired, do not hesitate to contact BTI Technical Service.

### INSPECTION, CLEANING AND DISINFECTION

To ensure a correct functioning of the oven it is recommended to inspect and clean the equipment regularly.

The cleaning frequency will depend on the use of the equipment and must be established by the user.

Whenever there is an accidental spill of plasma affecting the oven (inside, tray, outer case, front panel...) it is recommended to clean and disinfect the unit.

 **DANGER**  
*Before cleaning the oven, turn OFF the equipment and disconnect it from the power supply.  
 If the unit has been used recently, wait until it is cold to proceed with the cleaning process.*

### CLEANING

- When cleaning the tray, extract it from the unit and proceed with the cleaning protocol.
- Use a NEUTRAL DETERGENT for cleaning. DO NOT USE ALKALINE DETERGENTS since they will damage the tray.
- Clean the unit with a NEUTRAL DETERGENT dampened cloth.
- Rub gently until it is clean, rinse the cloth with water and then wipe again to remove all the detergent.
- Dry with an absorbent cloth.
- Do not use abrasive products.

### DISINFECTION

- If disinfection is required, make sure first to clean properly and rinse it well with water to eliminate any type of dirt and then proceed to disinfect the unit by using a 70% Ethanol solution.
- Rub gently all the affected surfaces with a 70% Ethanol solution dampened cloth.
- Do not dry the unit with a cloth; wait until the ethanol solution evaporates by itself.

 **WARNING**  
*It is necessary to check regularly the state of the unit and tray to clean them whenever required and to avoid accumulating substances that may affect the functioning of the unit.*

 *Follow the recommended precautions of the cleaning and disinfection products you are going to use.*

 *Do not sterilize the tray by autoclave; the sterilization process will harm the tray.*

### CHANGING FUSES

The fuses are located in the base of power connector. To replace them, first unplug the power cable and use a screwdriver to lever the fuse box cover. Extract and change the deteriorated fuses for new ones of the same amperage, which will be indicated on the characteristics label of the unit (picture 3) and in section "8. Technical Characteristics" of this manual.

 **DANGER**  
*Before changing the fuses, turn OFF the equipment and disconnect it from the power supply.  
 If the unit has been used recently, wait until it is cold to proceed with the changing process.*

### TECHNICAL ASSISTANCE

Electrical maintenance or repairs operations must always be carried out by specialised personnel authorised by BTI Technical Service.

If the unit has been used with pathogenic materials and requires technical assistance, the user must complete the procedure for cleaning and disinfection according to section 6.1 of this manual, to ensure the safety of the Technical Service staff before performing any repair or maintenance work.

- *In accordance with the current legislation EN 61010-2-010 Sec. 5.4.5, we recommend you contact BTI Technical Service for the opportune maintenance and verification work on the unit with a recommended annual frequency.*
- *This verification activity will be carried out by specialised technical staff authorised by BTI for the Maintenance and Technical support service.*
- *Remember that the correct functioning of the unit largely depends on the maintenance work and periodical verifications established by BTI.*

## 7. TROUBLESHOOTING

- *The unit is connected to the power supply, the master switch is ON and the status indicator does not light up.*

Change the fuses according to section 6.2 of this Manual.

- The unit is connected to the power supply, the master switch is ON and both the status and program indicators are blinking.

Contact BTI Technical Service. The unit may have a fault.

### 8. TECHNICAL CHARACTERISTICS\*

- Power supply (\*\*): 100 - 240V
- Maximum frequency: 50/60 Hz
- Maximum current assigned: 1.5A
- Fuse (type T): 2 x T 2 AL (5 x 20mm)
- Class (according to CISPR 11) (\*\*\*) (\*\*\*): 1B
- Working temperature: 37±3°C or 56±3°C
- Dimensions: (W/D/H): 241 / 221 / 245.5 mm
- Weight: 5.6 Kg.
- Warranty: 24 (months)
- Pollution grade: I

- Recommended conditions of use:
- Temperature margin: 10°C - 37°C
  - Relative humidity: <75%

- Transport and storage conditions:
- Humidity: < 90%
  - Temperature: -18°C a 40°C

- \* These technical characteristics are for:
- Ambient temperature between: 10°C and 37°C.
  - Grid margin: ± 10%.
  - Relative humidity: <75%
  - For installation category: I.
  - Pollution class: I
  - Indoor use
  - Altitude: up to 2000m

\*\*The grid must be earthed and the unit must be connected to it.  
 \*\*\*It is authorised for use in domestic establishments, like hospital or clinics under the responsibility of a medical professional.

### 9. SYMBOLS

Information present on the box label:

 DANGER	Caution, risk of danger Indicates to the operator that this method of use may cause physical harm.	 WARNING	Caution, see attached documentation. Indicates to the operator a method of use that could damage the unit.
	Biohazard.		Attention: Hot surface.
	Danger, risk of electric shock.		Conductor terminal for protection.
	CE Marking according to Council Directive 93/42/EEC corrected by the new Directive 2007/47/EC concerning medical devices (European requirement)		Selective collection of electrical and electronic devices. The equipment must be obligatorily disposed of separate from domestic waste.
	Serial number		Item number
	Batch number		Manufacturer
	See instructions for use		US federal law restricts this device to sale by or on the order of a licensed healthcare practitioner.

### 10. UNIT DISPOSAL

 The presence of this symbol on the product, its accessories or the packaging indicates that the product must not be treated as domestic waste. Dispose of the unit at your nearest collection point for recycling of electrical and electronic waste. In the European Union and other European countries there are different collection systems for used electrical and electronic products. Ensuring this product is disposed of correctly will help avoid potential risks for the environment and the health of people that could occur if the product is not disposed of properly. Recycling of materials helps preserve natural resources. Therefore, do not dispose of used electrical and electro-

nic equipment or batteries in the domestic waste. For more detailed information on the recycling of this product, contact your municipal office, the domestic waste collection service or the establishment where you acquired the product.

## 11. SAFETY INFORMATION. WARNINGS

### WARNINGS

In order to maintain the operator's safety and ensure the correct functioning of the unit, follow carefully the notes below as well as all the warnings in this manual.

 **DANGER:** For greater safety, the unit is designed to partially prevent spills from accessing to electrical parts. Ensure the containers are not full of plasma to the limit of their capacity to avoid potential spills when placing them in the unit and be sure to close them properly.  
Do not lean on the unit or place objects on it while it is functioning.  
If the unit is used in a manner other than that specified by BTI, user safety and unit protection may be compromised.

### SAFETY

This unit, like any other electrical equipment, is a potential source of electromagnetic radiation. It meets the requirements for electromagnetic compatibility, emission and immunity for electrical equipment for measurement, control and use in laboratories. (IEC 61326-1).

In any event, you must evaluate the electromagnetic environment before using the unit.

The unit includes safety systems for preventing overheating in the event of a fault.

 **WARNING** Do not use the unit near sources of strong electromagnetic radiation (for example, unshielded intentional RF sources), as they may interfere with the proper operation of the unit.  
Check the unit's power supply requirements, so that it is not connected to an unsuitable source (check section 8. Technical Characteristics).  
For better conservation of the unit, carefully follow the instructions in this manual.

 **DANGER, HOT SURFACE**  
*The equipment is designed to heat and maintain at 37°C or 56°C Plasma Rich in Growth Factors (PRGF), within the scope of Endoret® technology. Therefore there are surfaces exposed to the user that may be hot during its use.*

## 12. WARRANTY

BTI guarantees this unit against all manufacturing defects for a period of 2 years, from the purchase date, in the following cases:

- ✓ Covering any manufacturing defect, including the labour necessary to locate and replace the defective parts in our Technical Service.
- ✓ This guarantee DOES NOT COVER faults that the BTI Technical Service deems to have been caused by incorrect installation, mistreatment, improper use or manipulation by people other than the BTI Technical Service.
- ✓ Any acknowledgement by BTI of compensation due to direct or indirect harm of any nature suffered by people or objects is expressly excluded.
- ✓ Any unit that has had the characteristics label removed or altered is not covered by this warranty.

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