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Stand der Information: 2017-05

Viscoelastic solution
for ocular surgery

INSTRUCTIONS FOR USE

VISIOL®

Sodium hyaluronate 2.0%. Viscoelastic solution for ocular surgery. Sterile by moist heat.

Composition:

1 ml isotonic solution (pH 6.8–7.5; osmolality: 270–330 mOsmol/kg) contains 20.0 mg sodium hyaluronate from fermentation and sodium chloride, disodium phosphate, sodium dihydrogen phosphate, mannitol and water for injections.

Indications:

VISIOL® is indicated for use as a surgical aid in the surgery of the anterior segment of the eye such as cataract surgery, intraocular lens (IOL) implantation, corneal transplantation surgery and glaucoma filtering surgery.

Contra-indications:

VISIOL® should not be used in patients with ascertained hypersensitivity to any of its constituents.

Interactions:

Avoid using VISIOL® with instruments sterilised with quarternary ammonium salt solutions.

Undesirable effects:

Sodium hyaluronate is a natural component of the extracellular matrix and is very well tolerated in the eye. To date, there have been no reports of any postoperative inflammatory reactions that are causally related to the use of VISIOL®. A transient post-operative increase in intraocular pressure (IOP) has been reported with other sodium hyaluronate products when the product was only partially removed from the anterior chamber. Clinical trials have shown that VISIOL® did not cause clinically significant increase in IOP if some product remained in situ after the surgery. However, a transient slight increase in IOP after the use of VISIOL® cannot be completely excluded.

Dosage and administration:

The content and outer surface of the VISIOL® pre-filled syringe are sterile as long as the sterile pack remains unbroken. Take the pre-filled syringe out of the sterile pack, unscrew the Luer-Lock cap, attach a suitable needle (for example ophthalmologic needle 25 G) and secure by turning slightly. Remove any air bubble, if present, before injection.

Cataract surgery and IOL implantation: VISIOL® can be used at any stage of cataract surgery to create a deep anterior chamber, protect the tissues and facilitate the IOL implantation. Extrude the required amount of VISIOL® slowly and carefully into the anterior chamber through a cannula. VISIOL® may also be used to coat the surgical instruments and IOL before insertion. Additional VISIOL® can be injected during surgery, if needed.

Corneal transplant surgery: Remove the corneal button and fill the anterior chamber with VISIOL® until it is level with the surface of the cornea. Place the donor graft on top of VISIOL® and suture into place. Additional VISIOL® can be injected during surgery, if needed.

Glaucoma filtration surgery: Inject the required amount of VISIOL® slowly and carefully into the anterior chamber through a paracentesis when, performing the trabeculectomy. Additional VISIOL® can be injected during surgery, if needed.

Hornhauttransplantation: After removing the corneal button and filling the anterior chamber with VISIOL®, suture the corneal button onto the recipient cornea. Place the donor graft on top of VISIOL® and suture into place. Additional VISIOL® can be injected during surgery, if needed.

Glaucomfiltration: During trabeculectomy through a paracentesis, fill the anterior chamber with VISIOL® until it is level with the surface of the cornea. Place the graft on top of VISIOL® and suture into place. Additional VISIOL® can be injected during surgery, if needed.

Precautions:

The normal precautions associated with eye surgery should be observed to avoid intra- and/or post-operative increase in IOP: Do not overfill the ocular chamber with VISIOL®. The product must be administered with care and under close monitoring, particularly in patients with pre-existing glaucoma and in cases of glaucoma surgery. If intra-ocular pressure rises above normal after surgery, appropriate treatment should be provided. The product should be removed by irrigation and/or aspiration at the end of the procedure. Do not use if the pre-filled syringe or sterile pack are damaged. Any solution not used immediately after opening must be discarded. Otherwise the sterility is no longer guaranteed and this may be associated with a risk of infection. Store between 2 °C and 25 °C! Do not use after the expiry date indicated on the box. Please dispose of the used syringe and the packaging respecting the environment and in accordance with the local waste disposal regulations.

Characteristics and mode of action:

Sodium hyaluronate, the active principle in VISIOL®, is a polysaccharide which consists of repeating sequences of glucuronic acid and N-acetylglucosamine. It is present in the extracellular matrix, in particular in the vitreous humour, the synovial fluid and the umbilical cord. The sodium hyaluronate in VISIOL® is obtained by fermentation and is highly purified. It has an average molecular weight of 1.7 – 1.8 million daltons. A graphic representation according to ISO 15798 shows the rheological profile of VISIOL®. VISIOL® exhibits a pseudoplastic flow behaviour, i.e. the viscosity decreases when the shear rate is increased. The extrapolated zero-shear viscosity is approximately 60,000 mPas.

VISIOL® helps create and maintain anterior chamber depth and visibility at all stages of the anterior segment surgery and minimises interaction between tissues during surgical manipulation. VISIOL® protects intraocular tissues, such as the corneal endothelium, from damage due to the use of surgical instruments during surgery. VISIOL® may also be used to coat the surgical instruments and IOL before insertion. VISIOL® also protects corneal endothelium against damage caused by free radicals. VISIOL® preserves tissue integrity and provides good visibility during eye surgery. The mannitol is present in VISIOL® to help to maintain the rheological properties of the hyaluronic acid in VISIOL® during phacoemulsification in cataract surgery.

Presentation:

One pre-filled syringe of 20 mg/1.0 ml VISIOL® in a sterile pack.

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