
# ENDORET<sup>®</sup> OPHTHALMOLOGY KIT

<b>en</b>	.....	3
<b>es</b>	.....	9
<b>de</b>	.....	16
<b>fr</b>	.....	23
<b>it</b>	.....	30
<b>pt</b>	.....	37
<b>cs</b>	.....	44
<b>sk</b>	.....	50
<b>sv</b>	.....	56
<b>no</b>	.....	62

# ENDORET<sup>®</sup> OPHTHALMOLOGY KIT

## INTENDED USE

ENDORET<sup>®</sup> OPTHALMOLOGY KIT is a system for obtaining plasma rich in growth factors using autologous technology. Its formulation of eye drops can be used to enhance the regeneration of ocular tissue.

## PRODUCT DESCRIPTION

ENDORET<sup>®</sup> OPTHALMOLOGY KIT (BTI Ref. KMU11) contains disposable kits (see amount on label). Each one of the disposable kits contains the following trays.

### Tray 1: Blood Extraction Tray (BBT)

REFERENCE	PRODUCT DESCRIPTION	STERILITY
TB9	COLLECTION TUBES SODIUM CITRATE 3.8%. (8 units)	Sterile*
450085V1	SAFETY Blood Collection Set + Holder (1 unit)	Sterile

### Tray 2: Fractionation Tray (OFT)

REFERENCE	PRODUCT DESCRIPTION	STERILITY
PTD2	BTI PLASMA TRANSFER DEVICE (1 unit)	Sterile
TP10	BTI 10 ML PLASMA TUBE (2 unit)	Sterile

### Tray 3: Ophthalmic Application Tray (OAT)

REFERENCE	PRODUCT DESCRIPTION	STERILITY
GP31830	DROPPERS OF 1 ML (32 units)	Sterile
GP32111	LUER-TP10 CONNECTOR (2 unit)	Sterile
GP31817	21G DISPENSER NEEDLES (2 unit)	Sterile
GP31802	ACTIVATOR CAPSULE (1 unit)	Sterile
GP31803	ACTUATOR WHEEL (1 unit)	Sterile
GP31840	DROPPER RACK (1 unit)	Sterile
GP31841	DROPPER RACK COVER (1 unit)	Sterile

### Tray 4: Ophthalmology prevention Tray (OST)

REFERENCE	PRODUCT DESCRIPTION	STERILITY
TC	Conservation container (1 unit)	Non-sterile
BOLSA OF-TALMO	BAG OPTHALMOLOGY KIT (1 unit)	Non-sterile
324892	BD Micro-Fine™ sterile insulin syringe (1 unit)	Sterile
ACT1	ENDORET ACTIVATOR 1mL (1 unit)	Sterile*
8300013707	Henke-Ject® needle 14Gx1 ½'' 2,1x40mm (8 units)	Sterile
SLGP033RB	Millex-GP Syringe Filter Unit, 0.22 µm, polyethersulfone, 33 mm, gamma sterilized (2 unit)	Sterile
ET.16X55 BLANCA**	Patient labels. 5 units	Non-sterile

\*Medical device with non-sterile exterior. It is advisable to clean the exterior of the vial with sterile gauze with alcohol at 70%.

\*\* It is not a medical device.

This reference KMU11 also includes the instructions for use: MA130 and MA142 (as many units as disposable kits).

## CONTRAINDICATIONS:

- Do not inject PRGF-ENDORET<sup>®</sup> or ENDORET ACTIVATOR 1mL intravascularly.
- Do not use ENDORET ACTIVATOR 1mL in cases of hypersensitivity to calcium chloride.
- Do not use PRGF-ENDORET<sup>®</sup> in cases of hypersensitivity to calcium citrate.

## ADVERSE EFFECTS:

PRGF-ENDORET<sup>®</sup> related products do not have adverse effects when used as it is detailed in section: Intended use. Nevertheless, some general recommendations regarding PRGF-ENDORET<sup>®</sup> should be taken into account.

- Do not mix PRGF-ENDORET<sup>®</sup> with other products, except biomaterials that are mixed with the plasma to form grafts.
- Do not use PRGF-ENDORET<sup>®</sup> in patients with active infectious processes.
- Do not use PRGF-ENDORET<sup>®</sup> if the patient presents blood disorders related to platelet or coagulation alterations.

- Do not use PRGF- ENDORET® to treat patients with an infection, skin disease or tumour in the injection area.
- The contraindication related to the application technique of PRGF-ENDORET®, clinical or surgical, must be observed/taken into account by the clinician.
- A light diet is recommended during the days prior to the treatment with PRGF- ENDORET®.

## CAUTIONS

- Handle all biological samples and sharp objects for blood extraction (lancets, needles, Luer adapters and blood collection equipment) according to health and safety policies and installation procedures.
- Viruses (such as hepatitis B and C, HIV (AIDS) or other infectious diseases) can be transmitted by contact with biological samples. In the event of exposure to biological samples, seek appropriate medical care.
- Do not use the tubes if foreign material is detected inside.
- The reuse and/or sterilisation of these products may cause risk of infection and/or cross contamination.
- Always follow the blood extraction protocol of the centre.
- Children, pregnant and lactating women: no extra precautions are considered, however, the individual situation of each patient must be taken into account by the clinician before applying the PRGF-ENDORET® technique.
- The clotting times in the preparation of PRGF-ENDORET® can be affected (increased) in patients treated with anticoagulants.
- The health/clinical professional is responsible for informing the patient of the contraindications, precautions and adverse effects related to the application technique of PRGF-ENDORET®.
- Any other additional precaution related to the application technique of PRGF-ENDORET®, clinical or surgical, must be observed/taken into account by the clinician.
- Contact your distributor if the package is opened or altered.

NOTE: The devices should be accepted only if the factory packaging and labeling arrive intact.

NOTE: If hemolysis or coagulated blood is observed after the centrifugation, discard the tubes.

## ASEPTIC PRACTICE

- Follow the aseptic operation protocols of your centre to minimise the possibility of contaminating the plasma fractions obtained through the PRGF-ENDORET® technique.
- Carry out the fractionation and activation of the plasma in a clean environment or in the surgery.
- Wear gloves during blood extraction and while handling tubes to minimise the risk of exposure.

## WASTE DISPOSAL

- Observe the general hygiene directives and legal regulations governing the proper disposal of infectious material.
- Wear gloves to reduce the risk of infection.
- Full collection tubes or medical devices containing contaminated blood must be collected in containers suitable for potentially infectious waste.
- Waste is usually disposed of either by incineration or by autoclaving

(steam sterilisation).

## KIT STORAGE PRECAUTIONS:

- Storage temperature: store between 4 and 25 °C.
- Avoid direct exposure to sunlight.
- Keep in a dry atmosphere.
- Do not use after the expiry date. Check the expiry date on the label of the single-use kit.

 **WARNING:** Exceeding maximum storage temperature may cause a deterioration in tube and vial quality (i.e. loss of vacuum, evaporation of liquid additives, colour change, etc.).

 **WARNING:** Exposure to humidity can affect the integrity of the tray.

## INSTRUCTIONS FOR USE:

NOTE: BTI regularly carries out training courses for the correct application of the ENDORET® technique.

### 1. EXTRACTION OF BLOOD

-  Gloves must be worn throughout the entire process.
-  Blood should be used within 4 hours of extraction. Do not extract blood if you plan to use it more than 4 hours after extraction.

1. Open the Blood Extraction Tray (BBT) and take out the material contained therein.
2. Identify COLLECTION TUBES SODIUM CITRATE 3.8% (TB9) for each patient.
3. Apply a compressor (max. 1 minute). Prepare the venipuncture site with an antiseptic solution. Do not touch the venipuncture area after cleaning.
4. Extract the SAFETY Blood Collection Set + Holder (450085V1) from its sterilised package.
5. Remove the cap from the needle.
6. Perform the venipuncture with the patient's arm in a downward position.
7. Push the tube against the holder and the valve of the needle, piercing the rubber diaphragm. Centre the tubes properly in the stand when inserting the cap to prevent lateral perforation and subsequent premature vacuum loss.
8. Remove the compressor as soon as the blood appears in the tube. Do not allow the contents of the tube to touch the cap during the procedure. Do not invert the tube and always keep it in an upwards position squeezing it with your thumb to ensure complete vacuum aspiration. The vacuum in the tube enables the flow of blood to stop automatically.
9. When the first tube is full and the blood flow stops, remove it gently from the holder.

N.B.: If no blood flows into the tube or if the flow stops before an adequate sample is obtained, the following steps are suggested to complete satisfactory blood extraction.

- a. Push the tube forward to pierce the cap. Keep it in position at all times, squeezing with your thumb to ensure complete vacuum extraction.
- b. Confirm the correct position of the needle in the vein.
- c. If the blood still fails to flow, remove the tube and place a new one in the holder.
- d. If this action still fails to give satisfactory results, remove the needle and discard it. Repeat the procedure from step 1.

10. Gently invert the tubes 4 to 6 times immediately after blood extraction to achieve a suitable mixture of additive and blood. Turn the full tube upside-down and back again. This is a complete inversion.

**!** Do not shake the tubes. Energetic mixing can cause foam, hemolysis and activation of the platelets. Insufficient mixing of the blood/additive can also produce undesirable results for the ENDORET® technique. Ensure the blood has been mixed with the additive before continuing.

11. Repeat steps 7 to 10 with the rest of the tubes required.

12. As soon as the blood stops flowing into the last tube, remove the needle from the vein, applying pressure to the puncture area with a dry sterile swab until bleeding stops.

**!** After venipuncture, the top of the cap may contain traces of blood. It is advisable to clean the septum of TB9 with sterile gauze with alcohol at 70%.

13. Discard the used needle with the holder in accordance with section WASTE DISPOSAL.

## 2. CENTRIFUGATION:


**!** Once full, the tubes must be centrifuged immediately after the blood extraction or within a maximum of one hour. The blood samples must not be refrigerated for any reason.

**!** TB9 SHOULD ONLY BE CENTRIFUGED IN A BTI SYSTEM MODEL CENTRIFUGE. In this way, we ensure that the plasma fractions obtained by the ENDORET® technique are correctly separated. Follow the specific instructions of use for the BTI System centrifuge.

**After centrifugation, the blood is separated into three components:**

- Plasma rich in growth factors: this is the column of plasma (yellowish) that contains most of the platelets, distributed according to a gradient of increasing concentration (this means that the number of platelets is lower at the top of the tube and increases in a downward direction).
- White cells or leukocytes: this is the fine whitish layer deposited just above the red blood cells (also called buffy coat).
- Red series, erythrocytes or red blood cells: this is the red column occupying the bottom of the tube.

Figure 1:

**!** If hemolysis (redness in the plasma) is observed after the centrifugation, discard the tubes.

**!** If centrifugation is not performed immediately after extraction, it is advisable to repeat step 1.10 (invert the tubes gently) before

centrifugation.

**!** When the tubes are placed on the stand after centrifugation, it is essential to keep the tubes vertical when handling them and refrain from making brusque movements of shaking them to prevent the different components being mixed.

## 3. GENERATION OF VACUUM IN THE BTI TUBE (TP10):

**!** Gloves must be worn throughout the entire process.

1. Open the Fractionation Tray (OFT) and take out the material contained therein.

N.B.: It is advisable to label the TP10 with the patient labels which are in the Ophthalmology prevention Tray (OST). Do so longitudinally so the scale is not hidden.

2. Move the plunger of the TP10 back to create the correct vacuum inside the tube, avoiding any sudden movements or the rotation of the axis (a). To ensure the vacuum is generated, hold the clamping ring firmly next to the lateral flanges preventing it from coming loose (a and b).

**!** If the clamping ring becomes loose, reinsert it in the tube with the tab until you hear a "click".

3. When the plunger is as far back as it goes (you must hear a "click"), turn the plunger anticlockwise to separate it from the tube.

4. To produce the vacuum in the second TP10, remove the clamping ring next to the lateral flanges of the TP10, insert the plunger until it fits inside the vacuum plunger, rotate it clockwise at least ¼ turn until its position enables you to reinsert the clamping ring. Insert the clamping ring so that the flange faces up until you hear a "click" (b). Obtain the vacuum in the TP10 as described previously and remove the plunger from the tube.

**!** To ensure the correct vacuum generated, use the tube immediately.

## 4. FRACTIONATION USING THE PLASMA TRANSFER DEVICE (PTD2)

**!** Gloves must be worn throughout the entire process.

1. Take the needles Henke-Ject needle 14G (8300013707) from the OST tray and insert them into the centre of the septum of each of the centrifuged TB9. Take care not to mix the red cells with the plasma.

2. Remove the PTD2 from the OFT tray carefully.

3. Remove the black blocking system placed on one side of the PTD2

by lifting the ends until they touch each other and then removing it.

**!** **USEFUL LIFE IN STORAGE:** Once the blocking system is removed, the device is operational for one hour.

4. Remove the protective tube from the cannula.
5. Insert the TP10 in the back of the PTD2 guiding the fins of the tube in the grooves at the back of the PTD2. Once inside, make a quarter turn clockwise to adjust its position.

**!** **DO NOT ASPIRATE AIR.** Pressing the fuchsia button while the tip of the PTD2 is in the air will cause the TP10 to lose vacuum.

**!** **REMOVE THE BLOCKING SYSTEM.** If you place the TP10 in the PTD2 without removing the blocking system, the vacuum in the tube will be lost.

**!** **DO NOT REMOVE THE CAP OF THE TP10.** Keep the cap on the TP10 during the whole fractionation process to ensure the correct vacuum is maintained.

6. Insert the cannula of the PTD2 through the Henke-Ject needle 14G (8300013707) to the most superficial area of the plasma fraction. Keep the tip of the cannula of the PTD2 just beneath the surface of the plasma (around 1-2 millimetres) and in contact with the inner wall of the tube.
7. Gently press the fuchsia button to begin the extraction of the plasma. During the process, accompany suction with a vertical downward movement to keep the cannula just below the surface of the fraction you wish to extract at all times. Aspirate slowly to avoid turbulence.

**!** **PRESS THE BUTTON GENTLY:** if you press the fuchsia button too hard at the start of the aspiration, this will make it impossible to control the aspiration speed, which can cause turbulence.

8. Aspirate all the plasma fraction (Vt) contained in the TB9. Do not aspirate the erythrocytes, nor the leukocytes (buffy coat).
9. Repeat steps 6 to 8 for each TB9 with the same PTD2 and TP10 until it is completely filled with plasma.
10. Once the tube has been filled, remove the safety blocking ring while the tip of the PTD2 is still inside the TB9 and remove the TP10 from the PTD2 with a slight anticlockwise turn.

**!** **REMOVAL OF THE CLAMPING RING:** If you do not remove the clamping ring before removing the TB9 from the PTD2, air could get into the TP10.

11. Insert the second TP10, in which you have previously created the vacuum, in the PTD2.
12. Repeat steps 5 to 9 until the second TP10 is full of plasma.
13. Discard the single-use PLASMA TRANSFER DEVICE (PTD2). It is NON-REUSABLE.

#### 5. PLATELET ACTIVATION PROTOCOL

**!** Gloves must be worn throughout the entire process.

N.B.: Once the plasma has been prepared, activate it within a maximum of 4 hours after blood extraction.

**!** Vial with non-sterile exterior. It is advisable to clean the exterior of the vial with sterile gauze with alcohol at 70%.

1. Activation of the plasma is carried out by adding a precise quantity of ENDORET ACTIVATOR 1mL (ACT1) and maintaining a temperature of 37 °C throughout the activation process. BTI recommends its

plasma heating unit. (PLASMATERM).

2. Gently turn the TP10 upside down three times and calculate the volume of plasma it contains using the scale graduated in mL (volume) on the outside of the tube.
3. Add the necessary volume of ACT1 using the Micro-Fine 0.5mL sterile insulin syringe (324892), according to the volume of plasma to be activated, at a rate of 0.02 mL (2 units on the scale of the Micro-Fine syringe) of ACT1 per millilitre of plasma.

N.B.: Bear in mind that the scale on the Micro-Fine syringe is marked in UNITS, not millilitres. Conversion is essential to avoid errors.

Volume of plasma (ml)	Volume of ACT1 (mL)	Volume of ACT1 (Units marked on the Micro-Fine 0.5mL sterile insulin syringe)
1	0.02	2
2	0.04	4
4	0.08	8
6	0.12	12
8	0.16	16
10	0.20	20

**!** The user should respect the amounts of ACT1 to be added, otherwise the clotting times will be affected. Do not use more ACT1 than is recommended. The addition of more ACT1 will not hasten the clotting process.

4. Pierce the septum of TP10 with the syringe and add the necessary amount of activator. Gently invert the TP10 to mix the activator with the plasma.
5. Discard the syringe and the surplus ACT1.
6. When performing the regular treatment, keep the activated plasma at 37 °C for 40 minutes so that the coagulation and retraction of the plasma can occur. BTI recommends its PLASMATERM H OVEN.
7. If you wish to perform the immunosafe treatment, continue as follows in accordance with the instructions for use of the ENDORET® PLASMATERM H OVEN (MA074) and PLASMATERM H-PLUS OVEN (MA156):

- **If you have a PLASMATERM H and a PLASMATERM H-PLUS OVEN**
  - Keep the plasma in the PLASMATERM H oven for 40 minutes at 37 °C.
  - Then insert the TP10 in the PLASMATERM H-PLUS OVEN for 1 hour with programme 2 selected and stabilised.
  - Continue with step 6 of the protocol "6. Preparation of the eye drops".
- **If you only have a PLASMATERM H-PLUS OVEN**
  - After keeping the plasma in programme 1 for 40 min at 37 °C, without removing the TP10 from the PLASMATERM H-PLUS, press the button "PROG" once.
  - At this point, programme 3 will start to blink for 1 hour. After this time, the oven will notify that the treatment has finished.
  - Continue with step 6 of the protocol "6. Preparation of the eye drops".


N.B.: It is important to remove the TP10 from the PLASMATERM H-PLUS immediately after the treatment is finished.

## 6. PREPARATION OF THE EYE DROPS:


 Gloves must be worn throughout the entire process.

N.B.: Once the plasma has been activated, prepare the eye drops within a maximum of 4 hours after blood extraction.

1. Open Ophthalmic Application Tray (OAT). Take the ACTUATOR WHEEL (GP31803) with gloves and place it on the plunger turning it anticlockwise (c). Insert the set in the plunger of the TP10 by turning it clockwise a little (d). Place the plunger forming a "T" between the grooved plunger and the lateral flanges of the TP10 (e).

2. Adjust the wheel moving it up to the top of the tube (f) and place the ACTIVATOR CAPSULE (GP31802) (g). Remove the lower cap of the TP10 (h) and fit the LUER-TP10 CONNECTOR (GP32111) and press firmly (i). Finally, fit the Millex-GP Syringe Filter Unit, 0.22 µm, polyethersulfone, 33 mm, gamma sterilized (SLGP033RB) (j).

 A greater pressure of the wheel against the lateral flanges makes the "plunger-wheel" assembly come out of the TP10.

3. Fit the TP10 upright with the filter up and turn the wheel until the system is drained. To do this you must see that the plasma completely moistens the Millex-GP filter.

4. Fit the 21G DISPENSER NEEDLE (GP31817) in the Millex-GP filter in the conical Luer connection. Remove the cap from the needle.

5. Position the DROPPER RACK (GP31840) so that the lower part of the rack faces up and the microbiological control dropper is located farthest from you (k).

6. Insert the needle directly through the centre hole of the first DROPPER OF 1 mL (GP31830), located in the lower right part of the dropper rack (l). Turn the wheel until the plasma begins to come out. The wheel should be rotated 2 marks until the second mark occupies the same position as the original mark of the wheel (m), which is equivalent to 0.6 mL.

 Take care when turning the wheel up to the mark to avoid overdosing.

7. After turning the wheel, wait 4-5 seconds until the PRGF-ENDORET® stops flowing from the needle and it is no longer in contact with the PRGF-ENDORET® to be able to remove it from the dropper rack.

 Gently and gradually remove the dispensing needle from the dropper.

8. Repeat steps 6 and 7 for the first 16 droppers.

9. After completing the 16 first droppers, pierce the microbiological control dropper and fill with 0.6 mL (n).

10. Without removing the needle from the microbiological dropper, remove the activator capsule and the plunger-wheel set and dispose of the TP10 with the Luer connector, the Millex-GP Unit and the dispensing needle.

11. Repeat steps 1-7 until the remaining droppers of the rack are filled with the second TP10 of plasma.

12. Finally, apply the necessary volume of eye drops from the second TP10 in the dropper for microbiology to fill it, so this dropper contains PRGF-ENDORET® from both TP10. Finally, dispose of the TP10 with the Luer connector, the Millex-GP Unit and the needle, the capsule and the plunger.

13. Break the tab of the DROPPER RACK COVER (GP31841) and remove the dropper you are going to use for the microbiological control, if this is required.

## 7. HOW TO USE THE MILLEX-GP SYRINGE FILTER UNIT, 0.22 µm, POLYETHERSULFONE, 33 MM (SLGP033RB)

This section provides the warnings and precautions, as well as the steps that must be followed to use the Millex-GP Syringe Filter Unit, 0.22 µm, polyethersulfone, 33 mm (SLGP033RB).

### WARNINGS:

- To guarantee sterility, do not use this product if the packaging is damaged.
- Do not use this product as an in-line filter for the intravenous administration of liquids; it has not been designed for prolonged continuous use.
- Do not use with syringes smaller than 10 mL, because they could reach pressures higher than the maximum nominal pressure, which can damage the filter or cause injury to staff.

### CAUTIONS:

- Do not use the Millex® filtration unit to filter liquids at temperatures above 45°C.

- Do not use the Millex® filtration unit to filter emulsions or suspensions, because it is not designed for this purpose.
- Do not use the Millex® unit to filter solutions containing an amount equal to or less than 5 mg of pharmacologically active substances, except when bond studies have been carried out.
- Do not re-sterilise or reuse this filtration unit, as Merck Millipore Ltd. cannot guarantee its sterility, integrity or operation beyond single use.

**PROCEDURE FOR USE:**

See section 6 "Preparation of the eye drops".

**8. STORAGE BEFORE DELIVERY TO THE PATIENT:**

Once all of the droppers have been filled, the product must be delivered to the patient as soon as possible:

1. Complete all the required fields shown on the dropper rack label with the patient's data.
2. Open the BAG OPHTHALMOLOGY KIT (BOLSA OFTALMO) by tearing in the area just above the yellow band and insert the Conservation container (TC) and the dropper rack inside. To close the bag, roll up the open end of the bag and fold the yellow tabs that protrude from the protective bag.
3. Send the dropper rack with the droppers containing the preparation in a container that ensures that it is delivered at least at room temperature (below 20-25 °C) when shipment does not exceed 3 days or refrigerated (2-8 °C) when delivery is expected to be between 3 to 7 days. Relevant patient information (contraindications, side effects and precautions: MA142) shall also be provided.

N.B.: BTI recommends its COOLBOXes (CBOX, CBOXM), which maintain the temperature between 2-8 °C.

N.B.: the components needed to give the eye drops to the patient (labels, dispenser rack and bag) are included in the OST tray.

 After the preparation, the patient should store the eye drops in a freezer as soon as possible, at temperatures below -15°C.

**CLINICAL RESPONSIBILITY**

 For optimum preparation of the eye drops, the clinician must follow the instructions in this manual.

 Preserve the traceability of the preparation and the dispenser intended for microbiological control (if this control is required).

 Provide the patient with the Patient Information Manual MA-142 together with the preparation, in addition to indicating the preparation date, the recommended duration of the treatment and the dosage that must be followed.

 Advise the patient to strictly follow the instructions for use and the storage conditions described in the Patient Information Manual MA-142. Do not begin the treatment under any circumstances unless you have received confirmation from the clinician responsible.

**INFORMATION PRESENT ON THE LABEL:**

	Use by date	<b>LOT</b>	Batch code
	CE marking		Content sufficient for 3 uses
<b>REF</b>	Reference		Consult instructions for use
	Temperature limit		Keep out of sunlight
	Caution		Do not re-sterilise
	Do not re-use		Keep dry
	Do not use if package is damaged and consult instructions for use	<b>STERILE R</b>	Sterilised using irradiation
	Manufacturer	<b>STERILE EO</b>	Sterilised using ethylene oxide.
<b>Rx only</b>	US federal law states that this device may be sold under medical prescription only.	<b>STERILE I</b>	Sterilized using steam or dry heat
	Date of manufacturer	<b>MD</b>	Medical device
<b>UDI</b>	Unique device identifier		Non sterile
	Single patient multiple use		

NOTE: The devices should be accepted only if the factory packaging and labeling arrive intact. Contact your distributor if the package is opened or altered.

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